



*Leukoplast® -
the professionals'
choice*



With focused dedication in wound care and continued commitment to excellence and service Leukoplast® has benefitted our customers for over 100 years in Europe. Now in the US market we will launch the Leukoplast® brand, which will be marketed under the umbrella of the Hypafix® product line.

Leukoplast® Hypafix®
covering all fixation needs

Leukoplast®

Wound care in best hands



Trusted performance with Leukoplast® Hypafix®

When it comes to medical products, quality is the top priority. The Hypafix® range of wide-area fixation products is your trusted partner when reliable adhesion is required. Flexible in use and easy to apply.

The features:

- Fast fixation on different wound sizes and also on moving body areas
- Flexible adaptation for patient comfort
- Reliable adhesion
- Highly breathable materials to prevent skin maceration and to help the wound to heal
- Helps to prevent contamination
- Residual-free removal
- Radio-transparent
- Not made with natural rubber latex

Our top-seller in wide-area fixation


Leukoplast® Hypafix®

For wide-area fixation on frequently moved and highly contoured body parts

This highly conformable, very soft and breathable wide-area fixation was designed by our experts to comfortably secure wound pads and ointment dressings, as well as instruments, probes, catheters and measuring devices, especially on frequently moved and highly contoured parts of the body.

Stretchable: made for moving body parts

- Non-woven fabric
- Resistant to tearing and pulling



High wearing comfort

- Soft material
- Stretches with movement

Breathable

- Air- and moisture-permeable

Low-allergy

- Suitable for patients with sensitive skin
- Not made with natural rubber latex

Easy application

- Made for contoured body parts

Leukoplast® Hypafix® transparent

For waterproof wide-area fixation of wound dressings

This skin-friendly, ultra thin PU film was developed by our experts to protect reliably against water as well as bacteria and other contaminants.



Waterproof and almost invisible

- Ultra-thin PU film
- Protects reliably against water
- Allows bathing and hydrotherapy

Flexible adaptation to body contours

- Very thin, soft and flexible

Quick and easy application

- Intuitive application technique with distinctive red stripe

Excellent wearing comfort

- skin friendly yet strong adhesive




Supports infection prevention

- Protects against bacteria and contaminants

Breathable

Long-lasting adhesion**

- Better adhesion than competitors and lasts longer



** Source: 2014-08-14 ProDERM Study No 14.0083-23, BSN medical GmbH

Leukoplast® Hypafix® skin sensitive

For secure fixation and reliable adhesion on fragile or compromised skin

Leukoplast® has expanded its range to include Hypafix® skin sensitive, which is designed especially for use on patients with fragile or compromised skin. The soft silicone adhesive helps to protect against adhesive-related injuries.



Protects fragile skin

- Reliable adhesion with minimal impact on skin

High wearing comfort

Atraumatic removal

- easy to use and re-position

Extra gentle

- >94% of users confirmed their skin did not react at all*




Easy application

- Can be repositioned multiple times
- Can be separated if stuck together
- Residual-free removal

Skin-friendly silicone adhesive

Minimal impact on skin integrity



*Source: C1525, Clinical Study Report: Prospective, multicenter, observational study to evaluate the product characteristics of an adhesive fixation dressing. Fixomull gentle touch - A non-interventional study by Global Clinical Operations, BSN medical GmbH. (Fixomull is also marketed under Hypafix brand)

Product overview


Order now – your solution for reliable fixation



Hypafix®		
Size	BNR	UOM/Case
2 in × 11 yds / 5 cm × 10 m	7144301	24
4 in × 11 yds / 10 cm × 10 m	7144302	24
6 in × 11 yds / 15 cm × 10 m	7144303	24
2 in × 2 yds / 5 cm × 2 m	7193921	36
4 in × 2 yds / 10 cm × 2 m	7193919	12
6 in × 2 yds / 15 cm × 2 m	7193920	12



Hypafix® skin sensitive		
Size	BNR	UOM/Case
2 in × 5.5 yds / 5 cm × 5 m	7996601	24
4 in × 5.5 yds / 10 cm × 5 m	7996602	12
4 in × 2.1 yds / 10 cm × 2 m	7996604	12



Hypafix® transparent		
Size	BNR	UOM/Case
2 in × 11 yds / 5 cm × 10 m	7237800	12
4 in × 11 yds / 10 cm × 10 m	7237801	12
6 in × 11 yds / 15 cm × 10 m	7237802	12
4 in × 2 yds / 10 cm × 2 m	7237803	12

Please refer to the product label and / or package insert for full instructions on the safe use of these products.



Leukoplast®,
an Essity brand



BSN Medical Inc., an Essity company
5825 Carnegie Blvd., Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325
To order toll-free: Tel. (+1) 800 537 1063

63774 R1 ©2024 BSN Medical Inc. C24